

## Surgical Implants


USA and International Patent Protected

- For medium to large orbital floor defects
- Implant bends and holds desired shape without heating
- Implant provided with a sterile, non-porous transparent template made of the same Su-Por material
- Exclusive smooth template allows for easy repeated insertion and removal during fitting
- Available with a heat bonded barrier layer
- Outline and contour of the modified template may be used in modification and placement of implant
- Fixation with screws or sutures
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



### As easy as 1-2-3:

1. Shape Template and Rinse Clean
2. Transfer Shape to Implant
3. Use Template as Introducer



CAT#	Description	Thickness
4487	40mm Radius w Template	0.85mm
4488	40mm Radius w Template	1.5mm
4489	Membrane - 40mm Radius w Tempalte	1.0mm