



Surgical Implants

- Fits over the inferior orbital rim and extends superiorly and inferiorly medial to the inferior orbital nerve
- Un-matched shape retention and suturability
- Fully porous material for tissue ingrowth
- Fixation screws or sutures can be placed in any desired location
- Easy to trim with scissors, scalpel, or fluted burr
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4191	Left	25mm	26mm	2.5mm
4192	Right	25mm	26mm	2.5mm

