



- Surgeon can manipulate the height and projection of helical rim to match the patient's contralateral ear
- Provides structural base for temporal parietal fascia flap and skin grafts
- Vascular tissue flap required to prevent late exposure
- Implant bends and holds desired shape without heating
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



**Ear
Base**



**Helical
Rim**



CAT#	Description	A	B
4099	Ear Base - Left	32mm	63mm
4100	Ear Base - Right	32mm	63mm
4101	Helical Rim - Left	39mm	63mm
4102	Helical Rim - Right	39mm	63mm

