




Sellar Floor

Designed with Murali Guthikonda, M.D., F.A.C.S.

- Designed to repair the sellar floor following transsphenoidal surgery.
- Provides a bolster for the contents of Sella turcica preventing empty sella syndrome.
- Does not require fixation as surrounding tissue holds position until tissue integration occurs.
- Available in two sizes with the large having three small tabs oriented to facilitate handling and placement while the regular has a single tab.
- Retained Shape Memory™ provides malleability by hand.



CAT#	Description	A	B	C	D
4126	Regular	20mm	20mm	2.5mm	0.45mm
4129	Large + Membrane	40mm	40mm	2.5mm	0.73mm

