



Surgical Implants

Designed with Michael A. Burnstine, M.D.

- Acts to lift and fill the eyebrow
- “Minimally invasive brow suspension technique succeeds because of elevation of the ROOF fat with porous scaffold to allow tissue ingrowth to help stabilize, elevate, and fix the eyebrow” - Michael A. Burnstine, M.D.
- Composite microchannel-suture brow suspension implant
- Anchor and stabilize to the deep temporalis fascia at the hairline for temporal eyebrow elevation
- Easy to trim with scissors or scalpel
- MR Safe, Sterile, 5 Year Shelf Life



**Scan for Video
Demonstration**



CAT#	Description	A	B	C
4442	Small - Qty 1	10mm	5mm	0.45mm
4443	Large - Qty 1	21mm	9.5mm	0.45mm
4444	Small - Qty 2	10mm	5mm	0.45mm
4445	Large - Qty 2	21mm	9.5mm	0.45mm

